



Product Datasheet

Tnfsf9 Antibody (orb1671772)

Description	Tnfsf9 Antibody
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	Blocking, FC, WB
Immunogen	This antibody was raised by immunising LOU rats with soluble recombinant 4-1BBL:Fc fusion protein.
Target	Tnfsf9
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P41274
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed RAW...



Flow cytometry using the anti- 4-1BBL an...