



Product Datasheet

Havcr1 Antibody (orb1671784)

Description	Havcr1 Antibody
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, Functional Studies, IF, IHC
Immunogen	This antibody was raised by immunising female Lewis strain rats (Harlan Sprague-Dawley) subcutaneously with mouse TIM-1-Ig in complete Freund's adjuvant (CFA), followed by multiple 'boosting' with mouse TIM-1-Ig in PBS.
Target	Havcr1
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	Q5QNS5
Expiration Date	12 months from date of receipt.



Surface plasmon resonance (SPR) sensorgr...



ELISA of anti-TIM-1 antibody on TIM-1-Fc...