



## Product Datasheet

### Ly6g Antibody (orb1671788)

<b>Description</b>	Ly6g Antibody
<b>Reactivity</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	Depletion, FC, IHC, WB
<b>Immunogen</b>	This antibody was raised in Lou-M rats using a popliteal lymph node immunisation and fusion protocol as described by Martiniello et al (1994)
<b>Target</b>	Ly6g
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	batch dependent
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG kappa
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<a href="#">P35461</a>
<b>Expiration Date</b>	12 months from date of receipt.



Immunofluorescence staining of mouse spl...



Flow-cytometry using the Anti-Ly6G/Ly6C ...