


Product Datasheet

TIGIT Antibody [4D4], Rabbit IgG (orb1671807)

Description	TIGIT Antibody [4D4], Rabbit IgG
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, Functional Studies
Immunogen	Armenian hamsters were immunized with recombinant mouse TIGIT tetramers by a combination of s.c. and food pad immunization and booster injections. Draining lymph nodes were fused with Sp2/O-Ag14 cells, selected in HAT (hypoxanthine/aminopterin/thymidine) medium and supernatants were screened for specificity by ELISA and flow cytometry using TIGIT-transfectants.
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P86176
Expiration Date	12 months from date of receipt.



ELISA of anti-TIGIT antibody on TIGIT-Fc...