



Product Datasheet

Cd276 Antibody (orb1671831)

Description	Cd276 Antibody
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	Blocking, FC
Immunogen	This antibody was generated by immunizing Sprague Dawley rats with B7-H3-Ig, consisting of the extracellular domain (aa1-242 of mouse B7-H3) linked to the Fc portion of mouse IgG2a, emulsified in CFA (Difco). The antibody was selected by its reactivity to mouse B7-H3-transfected cells, but not to parental cells by flow cytometry.
Target	Cd276
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	Q8VE98
Expiration Date	12 months from date of receipt.



Binding curve of anti-B7-H3 antibody MJ1...



Immunofluorescence staining of fixed mou...