



Product Datasheet

Icosl Antibody (orb1671835)

Description	Icosl Antibody
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	Blocking, FC
Immunogen	This antibody was raised by immunizing SD rat with mouse B7h-transfected L cells and fusing immune splenocytes with P3U1 myeloma cells and screened for binding to mouse B7h-transfected NRK cells.
Target	Icosl
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	Q9JHJ8
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed mou...



Binding curve of anti-ICOSL antibody HK5...