



Product Datasheet

MHC I Antibody [R1-21.2], Rabbit IgG (orb1671839)

Description	MHC I Antibody [R1-21.2], Rabbit IgG
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC, IP, RIA
Immunogen	This antibody was raised by immunising Lewis 1F rats with mouse CH-1 lymphoma cells. Spleen cells from immunised rats were fused with X-63Ag8-653 myeloma cells to generate hybridomas.
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of mouse spl...



Western Blot using anti-MHC I antibody R...