



Product Datasheet

MHC II Antibody [P7/7], Rabbit IgG (orb1671843)

Description	MHC II Antibody [P7/7], Rabbit IgG
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, IP
Immunogen	This antibody was raised by immunising DA rats with PVG.R8 rat spleen cells. Splenocytes from immunised rats were then isolated, and fused to the myeloma line NSO/u to generate stable hybridomas.
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Expiration Date	12 months from date of receipt.



Flow-cytometry using the anti-MHC II ant...



Immunofluorescence staining of mouse spl...