

Product Datasheet

Pdpm Biosimilar Antibody (orb1671846)

Description	Pdpm Biosimilar Antibody
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	Blocking, ELISA, Epitope Tagging, IHC, WB
Immunogen	This antibody was raised by immunising rats with 14-residue synthetic peptide mpp3851, which corresponds to amino acids 38-51 (GDGMVPPGIEDKIT) of the platelet-aggregation-stimulating (PLAG) domain of mouse podoplanin. Spleen cells were then harvested, and fused to P3U1 cells to generate stable hybridomas.
Target	Pdpm
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG2a kappa
Clonality	Recombinant
Uniprot ID	Q62011
Expiration Date	12 months from date of receipt.