



Product Datasheet

PD-1 Antibody [RMP1-14], Rabbit IgG (orb1671853)

Description	PD-1 Antibody [RMP1-14], Rabbit IgG
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	Blocking, FC, IHC-Fr
Immunogen	The original rat hybridoma was generated by immunizing Sprague Dawley rats with mouse PD-1-transfected BHK cells. A P3U1 myeloma was used as the fusion partner.
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Expiration Date	12 months from date of receipt.



Flow-cytometry using anti-PD-1 antibody ...



Immunofluorescence staining of fixed mou...