



Product Datasheet

Fc gamma receptor III Antibody [3G8], Rabbit IgG (orb1671874)

Description	Fc gamma receptor III Antibody [3G8], Rabbit IgG
Reactivity	Human, Marmoset, Monkey, Primate
Conjugation	Unconjugated
Tested Applications	Blocking, Depletion, FC, IHC, IP
Immunogen	3G8 was prepared by immunizing mice with human PMN (polymorphonuclear leukocytes) cells.
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	Q9ULV2
Expiration Date	12 months from date of receipt.



Flow-cytometry using the anti-Fc gamma r...



Immunofluorescence staining of fixed U93...