


Product Datasheet

DDDDK-tag Antibody [M2.1], Mouse IgG1 (orb1671880)

Description	DDDDK-tag Antibody [M2.1], Mouse IgG1
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Interleukin 2 genetically fused to the peptided DYKDDDDK, corresponding to the DDDDK epitope tag.
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1 kappa
Clonality	Recombinant
Expiration Date	12 months from date of receipt.



Binding of Anti-DDDDK antibody to Multit...