

Product Datasheet

IL6R Biosimilar Antibody (orb1671882)

| | |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| Description | IL6R Biosimilar Antibody |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, FA |
| Immunogen | The original antibody was generated by CDR grafting of mouse PM-1 antibody onto human REI and NEW framework regions to create rhPM-1. |
| Target | IL6R |
| Preservatives | PBS with 0.02% Proclin 300. |
| Concentration | batch dependent |
| Storage | Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C. |
| Note | For research use only |
| Isotype | IgG1 kappa |
| Clonality | Recombinant |
| Uniprot ID | P08887 |
| Expiration Date | 12 months from date of receipt. |