


## Product Datasheet

### CD80 Biosimilar Antibody (orb1671883)

<b>Description</b>	CD80 Biosimilar Antibody
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	Blocking, ELISA, FC, IHC, IP, WB
<b>Immunogen</b>	Galiximab was prepared by immunizing cynomolgus monkeys with recombinant CD80 antigen. The variable regions of the light and heavy chains were then cloned by being incorporated into a cassette vector (N5LG1) containing human constant region genes and subsequently transfected into the Dg44 CHO cell line (Hariharan K, 2013)
<b>Target</b>	CD80
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	batch dependent
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG lambda
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<b>P33681</b>
<b>Expiration Date</b>	12 months from date of receipt.



Immunofluorescence staining of fixed  
Dau...