

## Product Datasheet

### CD22 Biosimilar Antibody (orb1671892)

<b>Description</b>	CD22 Biosimilar Antibody
<b>Reactivity</b>	Human, Primate
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	Blocking, FC, IP
<b>Immunogen</b>	This antibody was prepared by the humanization of LL2 (EPB-2), a murine anti-CD22 IgG2a raised against Raji Burkitt lymphoma cells. Murine sequences comprise 5-10% of the molecule, with the remainder being human framework sequences, which greatly reduces the potential for immunogenicity (Traczewski, 2010).
<b>Target</b>	CD22
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	batch dependent
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1 kappa
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<b>P20273</b>
<b>Expiration Date</b>	12 months from date of receipt.