


Product Datasheet

CCR4 Biosimilar Antibody (orb1671893)

Description	CCR4 Biosimilar Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Depletion
Immunogen	This antibody was prepared by humanization of a mouse anti-CCR4 mAb.
Target	CCR4
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P51679
Expiration Date	12 months from date of receipt.



Western Blot using anti-CCR4 antibody KW...