


Product Datasheet

C3 Biosimilar Antibody (orb1671894)

Description	C3 Biosimilar Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Blocking, FC, IHC
Immunogen	This antibody we generated by humanisation of the murine antibody R3, raised against EGFR purified from placenta membranes and A431 tumors.
Target	C3
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P00533
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed HeL...