

Product Datasheet

TNF Biosimilar Antibody (orb1671896)

Description	TNF Biosimilar Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Blocking, ELISA, IF, IHC
Immunogen	This antibody was generated by affinity maturation from the 2SD4 scFv, which in turn was selected on human TNF alpha by guided phage-display technology, using the murine anti-hTNF alpha antibody Mab32 as a template (c.f. also US Patent 6090382).
Target	TNF
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P01375
Expiration Date	12 months from date of receipt.