

Product Datasheet

VEGFA Biosimilar Antibody (orb1671901)

Description	VEGFA Biosimilar Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Blocking, ELISA, FC, IP, WB
Immunogen	Recombinant human VEGF.
Target	VEGFA
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1 kappa
Clonality	Recombinant
Uniprot ID	P15692
Expiration Date	12 months from date of receipt.