



Product Datasheet

Cd27 Antibody (orb1671910)

Description	Cd27 Antibody
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	Activation, FC, IF, IP
Immunogen	This antibody was produced by injecting Armenian hamsters intraperitoneally with a mixture of two murine CD27 tag-expressing Armenian hamster fibroblast clones (AR0524 and AR0530) in Freund's adjuvant
Target	Cd27
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P41272
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed mou...



Binding curve of anti-CD27 antibody LG.3...