


Product Datasheet

Phl p 1 Antibody [Clone 25], Rabbit IgG (orb1671914)

Description	Phl p 1 Antibody [Clone 25], Rabbit IgG
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	DOT, ELISA, WB
Immunogen	Recombinant Phl p 1, a grass pollen allergen.
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P43213
Expiration Date	12 months from date of receipt.



Binding curve of anti-Phl p 1 antibody C...