


## Product Datasheet

# Phl p 1 Antibody [Clone 25], Human IgE (orb1671915)

<b>Description</b>	Phl p 1 Antibody [Clone 25], Human IgE
<b>Reactivity</b>	Plant
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	DOT, ELISA, WB
<b>Immunogen</b>	Recombinant Phl p 1, a grass pollen allergen.
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	batch dependent
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Isotype</b>	IgE kappa
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<b>P43213</b>
<b>Expiration Date</b>	12 months from date of receipt.



Binding curve of anti-Phl p 1 antibody C...