


Product Datasheet

Bet v 1 Antibody [BV16], Rabbit IgG (orb1671918)

Description	Bet v 1 Antibody [BV16], Rabbit IgG
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	Crystallography, ELISA, WB
Immunogen	Bet v 1-A, a major allergen of the European white birch.
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P15494
Expiration Date	12 months from date of receipt.



Binding curve of anti-Bet v 1 antibody B...