


Product Datasheet

tau Antibody [Tau-5], Mouse IgG1 (orb1671925)

Description	tau Antibody [Tau-5], Mouse IgG1
Reactivity	Human, Mouse, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	EM, IF
Immunogen	Tau protein.
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1 kappa
Clonality	Recombinant
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed U25...