



Product Datasheet

GITR Antibody [YGITR765], Rabbit IgG (orb1671931)

Description	GITR Antibody [YGITR765], Rabbit IgG
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	Blocking, FC
Immunogen	nB2/6TG cells transfected with mouse GITR injected into DA rats.
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed mou...



Flow-cytometry using anti-GITR antibody ...