


Product Datasheet

CALM1P2 Antibody (orb1671934)

Description	CALM1P2 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Blocking, ELISA, FC, ICC, IF, IHC, WB
Immunogen	This antibody was produced by immunizing BALB/c mice with purified Hermes-1 antigen from KCA, a human B lymphoblastoid line that binds preferentially to mucosal HEV.
Target	CALM1P2
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P16070
Expiration Date	12 months from date of receipt.



Flow-cytometry using anti-CD44 antibody ...