


Product Datasheet

CALM1P2 Antibody (orb1671936)

Description	CALM1P2 Antibody
Reactivity	Bovine, Goat, Human, Mouse, Primate, Rabbit, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	Blocking, ELISA, FC, ICC, IF, IHC, IP, NeA, WB
Immunogen	This antibody was produced by immunizing a Fischer rat with human tonsil lymphocytes (10^7 cells intraperitoneally in Freund's complete adjuvant). A boost was given 2 weeks later (intravascularly in saline).
Target	CALM1P2
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P16070
Expiration Date	12 months from date of receipt.



Flow-cytometry using anti-CD44 antibody ...