



Product Datasheet

Sell Antibody (orb1671946)

Description	Sell Antibody
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	Blocking, FC, IHC, IP
Immunogen	C3H/eb mouse B lymphoma 38C-13.
Target	Sell
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P18337
Expiration Date	12 months from date of receipt.



Flow-cytometry using anti-CD62L antibody...



Flow-cytometry using the anti-CD62L (L-s...