

Product Datasheet

PPP1R1B antibody (orb167195)

Description Rabbit polyclonal antibody to PPP1R1B.

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ICC, IF, WB

Immunogen Recombinant fusion protein of human DARPP32(NP_115568.2).

Preservatives PBS with 0.01% thiomersal,50% glycerol,pH7.3.

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Isotype IgG


Clonality Polyclonal

MW 34KDa


Uniprot ID [Q9UD71](#)

Dilution Range WB: 1:500-2000, IHC-P: 1:50-200


Expiration Date 12 months from date of receipt.



Western blot analysis of extracts of var...



Immunohistochemical staining of rat lung...



Immunohistochemical staining of rat kidn...