


## Product Datasheet

### TNFSF4 Biosimilar Antibody (orb1671953)

<b>Description</b>	TNFSF4 Biosimilar Antibody
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	Blocking, FC
<b>Immunogen</b>	Oxelumab binds to human OX40L, the ligand of OX40, stimulatory receptor that is expressed on activated T cells, natural killer (NK) cells and natural killer T (NKT) cells.
<b>Target</b>	TNFSF4
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	batch dependent
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG kappa
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<b>P23510</b>
<b>Expiration Date</b>	12 months from date of receipt.



Flow-cytometry using the anti-OX40L anti...