


Product Datasheet

ITGAL Biosimilar Antibody (orb1671955)

Description	ITGAL Biosimilar Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Blocking, FC
Immunogen	Efalizumab binds to CD11a, a subunit of leukocyte function antigen-1 (LFA-1), which is expressed on all leukocytes.
Target	ITGAL
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P20701
Expiration Date	12 months from date of receipt.



Flow-cytometry using the anti-CD11a rese...