



Product Datasheet

C3 Biosimilar Antibody (orb1671957)

Description	C3 Biosimilar Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Blocking, ELISA, FC, Functional Studies, IF, WB
Immunogen	The parental mouse antibody was generated by immunizing BALB/c mice intraperitoneally with A431 cells in phosphate buffered saline (PBS). Later on the humanized version of the antibody was created by grafting CDRs of the murine antibody onto human constant regions.
Target	C3
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P00533
Expiration Date	12 months from date of receipt.



Flow-cytometry using the anti-EGFR antib...



Immunofluorescence staining of fixed A43...