



Product Datasheet

TNFRSF8 Antibody (orb1671982)

Description	TNFRSF8 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Activation, ELISA, FC, IF, IHC
Immunogen	Human CD30.
Target	TNFRSF8
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P28908
Expiration Date	12 months from date of receipt.



Flow-cytometry using the anti-CD30 (Tumo...



Immunohistochemical staining of human to...