



Product Datasheet

dsRNA Antibody [1D3], Rabbit IgG (orb1671985)

Description	dsRNA Antibody [1D3], Rabbit IgG
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	RDV-RNA-methylated bovine serum albumin complex.
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed cel...



Immunofluorescence staining of Zika viru...