



Product Datasheet

CD4 Biosimilar Antibody (orb1671990)

Description	CD4 Biosimilar Antibody
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	Blocking, ELISA, FC, IHC
Immunogen	Domain 1 of human CD4.
Target	CD4
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG4 lambda
Clonality	Recombinant
Uniprot ID	P01730
Expiration Date	12 months from date of receipt.



Western Blot using anti-CD4 antibody CE9...



Immunohistochemical staining of guinea p...