



Product Datasheet

ITGA5 Biosimilar Antibody (orb1671991)

Description	ITGA5 Biosimilar Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Blocking, ELISA, FC, IF, IHC
Immunogen	Heterodimeric $\alpha 5\beta 1$ Fc fusion protein
Target	ITGA5
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P08648
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed HeLa...



Flow-cytometry using the anti-alpha 5 be...