


Product Datasheet

CD40LG Biosimilar Antibody (orb1671995)

Description	CD40LG Biosimilar Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	Blocking, FC, WB
Immunogen	Human CD40L
Target	CD40LG
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P29965
Expiration Date	12 months from date of receipt.



Flow-cytometry using the anti-CD40L rese...