


Product Datasheet

TNNI3 Antibody (orb1672006)

Description	TNNI3 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Cardiac Troponin I peptide (KISASRKLQLKT).
Target	TNNI3
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1 lambda
Clonality	Recombinant
Uniprot ID	P19429
Expiration Date	12 months from date of receipt.



Western Blot using anti-Cardiac Troponin...



Immunofluorescence staining of fixed HeL...