



Product Datasheet

TUBB3 Antibody (orb1672015)


Description	TUBB3 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Human MBP (microtubule-binding protein).
Target	TUBB3
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG lambda
Clonality	Recombinant
Uniprot ID	Q13509
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed HeLa...



Flow-cytometry using the anti-Beta Tubul...



Western Blot using anti-Beta-Tubulin ant...