



Product Datasheet

TUBA1A Antibody (orb1672017)


Description	TUBA1A Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Human MBP (microtubule-binding protein).
Target	TUBA1A
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG lambda
Clonality	Recombinant
Uniprot ID	Q71U36
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human ce...



Western Blot using anti-Alpha-Tubulin an...



Flow-cytometry using the anti-Alpha Tubu...