



Product Datasheet

GITR Antibody [YGITR 860.103.5], Rabbit IgG (orb1672029)

Description	GITR Antibody [YGITR 860.103.5], Rabbit IgG
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Mouse GITR.
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Expiration Date	12 months from date of receipt.



Flow-cytometry using anti-GITR antibody ...



Immunofluorescence staining of fixed mou...