



Product Datasheet

MUC1 Antibody (orb1672037)

Description	MUC1 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	Human breast epithelial mucin.
Target	MUC1
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P15941
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of fixed MCF...



Western Blot using anti-breast epithelia...