



Product Datasheet

VEGFA Antibody (orb1672039)


Description	VEGFA Antibody
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	Blocking, ELISA, FC, IHC, IP, WB
Immunogen	Recombinant human VEGF.
Target	VEGFA
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1 kappa
Clonality	Recombinant
Uniprot ID	P15692
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat brai...



Immunohistochemical staining of rat skin...



Western Blot using anti-VEGF antibody.
R...