

Product Datasheet

RHO Antibody (orb1672041)

Description	RHO Antibody
Reactivity	Amphibian, Bovine, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, RIA, WB
Immunogen	Bleached bovine rod outer segment (ROS).
Target	RHO
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1 kappa
Clonality	Recombinant
Uniprot ID	P08100
Expiration Date	12 months from date of receipt.