



Product Datasheet

UBE2I Antibody (orb1672048)

Description	UBE2I Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide "GLFIIDGK", derived from the SUMO-conjugating enzyme UBC9, conjugated to KLH.
Target	UBE2I
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P63279
Expiration Date	12 months from date of receipt.



Flow cytometry using the Anti-UBE2I anti...



Immunofluorescence staining of Jurak cel...