



Product Datasheet

TAGLN Antibody (orb1672051)


Description	TAGLN Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide "AAEDYGVIK", derived from the actin cross-linking protein Transgelin, conjugated to KLH.
Target	TAGLN
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	Q01995
Expiration Date	12 months from date of receipt.



Flow cytometry using the antibody SAIC-3...



Immunofluorescence staining of HeLa cell...



Western Blot using anti-TAGLN antibody S...