



Product Datasheet

FSCN1 Antibody (orb1672052)

Description	FSCN1 Antibody
Reactivity	Human
Conjugation	Unconjugated
Immunogen	Peptide "TVEEAENIAVTSGVVR" derived from Fascin, which organises filamentous actin into bundles, conjugated to KLH.
Target	FSCN1
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	Q16658
Expiration Date	12 months from date of receipt.



Flow cytometry using the Anti-FSCN1 anti...



Immunofluorescence staining of HeLa cell...