


Product Datasheet

IL18 Antibody (orb1672062)

Description	IL18 Antibody
Reactivity	Human
Conjugation	Unconjugated
Immunogen	The KLH conjugated peptide "ISTLSCENK" derived from IL-18, an augmentor of natural killer cell activity in the spleen and stimulator of interferon gamma production in Th1 cells.
Target	IL18
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	Q14116
Expiration Date	12 months from date of receipt.



Western Blot using anti-IL-18 antibody S...



Flow-cytometry using the anti-IL-18 anti...