



Product Datasheet

BCL5 Antibody (orb1672065)


Description	BCL5 Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Peptide "VSLDVNHFAPDELTVK" derived from Heat shock protein beta-1, involved in stress resistance and actin organisation, conjugated to KLH.
Target	BCL5
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P04792
Expiration Date	12 months from date of receipt.



Western Blot using anti-HSPB1 antibody S...



Flow cytometry using the Anti-HSPB1 anti...



Immunofluorescence staining of HeLa cell...