



Product Datasheet

GSTP1 Antibody (orb1672067)

Description	GSTP1 Antibody
Reactivity	Human
Conjugation	Unconjugated
Immunogen	Peptide "YISLIYTNYEAGK" derived from Glutathione S-transferase P conjugated to KLH.
Target	GSTP1
Preservatives	PBS with 0.02% Proclin 300.
Concentration	batch dependent
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG kappa
Clonality	Recombinant
Uniprot ID	P09211
Expiration Date	12 months from date of receipt.



Flow cytometry using the Anti-GSTP1 anti...



Immunofluorescence staining of HeLa cell...